

File name: C:\Solid\Combi\Data\Kaufhold_shale\TS-12c_2022-03-16.144

SOLID REPORTING

Company name: ETH Zürich
Operator: MP

TEST INFORMATION

Test date: 16-03-22
Sample name: TS-12c
Comment: neues Material März2022
Sample mass (g): 7.011
Sample skeleton density (g/cm³): 2.8306
Test file name: C:\Solid\Instr2\Data\Kaufhold_shale\TS-12c_2022-03-16.P44
Combined with file (140): C:\Solid\Instr1\Data\Kaufhold_Shale\TS-12c_2022-03-16.P14
Mercury surface tension (N/m): 0.48
Mercury contact angle (°): 140.0
Test filling volume (mm³): 453.0 at P < 1 Pa
Starting hydr. press. of test (MPa): 0.0132
(Dil+Hg+Sample) weight (g): 147.004
Corrected weight (Dil.+Hg+Sample) (g): 147.031

ANALYTICAL CONDITIONS

Maximum test pressure (MPa): 400
Increase speed: 5 - 17 MPa/min
Increase method: Pascal Stepwise
Decrease speed: 3 - 11 MPa/min
Decrease method: Pascal Stepwise
Temperature of test (°C): 22.2
Mercury density @ test (g/cm³): 13.5409


PREPARATION

Wait time for filling (min): 10
Air pulses (N.): 1
Filling volume (mm³): 450.0
Preparation selection: AutoMatic

BLANK & DILATOMETER INFORMATION (Data from Pascal 440 file)

Blank date: 05-10-20
Blank filename: C:\Solid\Instr2\Blank\BK_D8080_2020-10-05_1058.P44
Comment: Blindwert 80_80
Blank Max pressure (MPa): 400
Blank Increase speed: 5 - 17 MPa/min
Blank Increase method: Pascal
Blank Decrease speed: 3 - 11 MPa/min
Blank Decrease method: Pascal
Dil. number: 8080
Dil. type: CD3
Dil. Cone length (mm): 30
Dil. Electrode gap (mm): 5
Dil. stem radius (mm): 1.5
Dil. weight (g): 59.865
Temperature of blank (°C): 22.8
Mercury density @ blank (g/cm³): 13.5394
Blank filling volume (mm³): 455 at P < 1 Pa
Starting hydr. blank press. (MPa): 0.013248
(Dil+Hg) weight (g): 176.068


SAMPLE RESULTS



RESULTS WITHOUT COMPRESSIBILITY CORRECTION

Total intruded volume (mm ³)	148.84		
Total intruded volume (mm ³ /g):	21.23	at pressure of MPa:	400.2106
Spec. Vol. by skeleton dens. (mm ³ /g) Vd:	26.56		
Bulk density (g/cm ³):	2.6326		
Envelope density (g/cm ³):	2.6328	at pressure of MPa:	0.014 Diam.(nm) 108147.5
Apparent density (g/cm ³):	2.7885	at pressure of MPa:	400.2106 Diam.(nm) 3.7
Porosity by skeleton density (%):	6.99	Calculated by skeleton density of	2.8306 g/cm ³
Porosity by Hg intrusion (%):	5.59		
Inaccessible porosity (%):	1.40		

TOTAL PORE SIZE



Derivative calculated with
Moving average points:
Smoothing factor:

11
0.95

Calc. press. range (MPa): 0.0 to 400.211
Total pore volume (mm³/g): 21.23
Total pore surface area (m²/g): 2.911
Average pore diameter (nm): 29.18
Median pore diameter (nm): 66.94
Modal pore diameter (nm): 61.77
Peak n. 2 Pore diamet. (nm): 53.19
Peak n. 3 Pore diamet. (nm): 54494.45
Span: 261.6378

with model :	Conical
Hydraulic (4V/A)	
at (mm ³ /g) :	10.62
at dV/dlogD (mm ³ /g, D) :	13.81
at dV/dlogD (mm ³ /g, D) :	13.62
at dV/dlogD (mm ³ /g, D) :	2.97
at Min. Vol. (%) 10.0	at Max. Vol. (%) 90.0

TOTAL PORE SIZE HISTOGRAM




TABLE TOTAL PORE SIZE

Pore diameter ranges (nm)	Medium point (nm)	Specific vol. (mm ³ /g)	Specific vol. (%)	Relative vol. (mm ³ /g)	Relative vol. (%)	Porosity (%)
120000.00-10000.00	65000.00	2.71	12.77	2.71	12.77	0.707
10000.00-100.00	5050.00	8.54	40.23	5.83	27.46	1.535
100.00-50.00	75.00	12.35	58.18	3.81	17.95	1.003
50.00-20.00	35.00	16.54	77.89	4.18	19.71	1.101
20.00-3.60	11.80	21.23	100.00	4.70	22.11	1.236